

Contec HydroPure SATWIPES System

6% hydrogen peroxide blended with purified water and Amplitude Theta polyester/cellulose wipes

Contec HydroPure is a blend of 6% hydrogen peroxide and purified water. Efficacious against bacteria, fungi, moulds, yeasts and spores, Contec HydroPure leaves little-to-no residue making it ideal for use in critical areas.

Contec HydroPure is provided as a point-of-use SATWIPES system. Each double bagged pack contains a canister of 100 x 54gsm Amplitude Theta polyester/cellulose dry wipes and a 500ml capped bottle of Contec HydroPure.

Simply add 500ml of HydroPure to the dry roll, 5 mins before use, and use the impregnated wipes within 7 days.

Contec HydroPure is 0.2 micron filtered, and filled in a Grade A environment. The sterile product is filtered into pre-irradiated components. The dry wipe canister is irradiated at no less than 25 kGy.



FEATURES	BENEFITS
Impregnated wipes at point-of-use	• Convenient and controlled application of disinfectant
Sporicidal activity	• Fast and easy to use
Hydrogen Peroxide breaks down to water & oxygen	• Very low residue so suitable for critical areas • Saves time since no residue removal required
100 polyester / cellulose wipes in a canister	• Ideal for a Grade A / B area • Convenient to use and allows controlled application
Filtered to 0.2 microns in a Grade B cleanroom	• Ensures the product is free from contamination and particulates
Provided sterile	• Suitable for Grade A and B cleanrooms
Not classed as corrosive	• Good operator acceptability as low odour and hazard
Double bagged in linear tear packaging	• Each bag is easy to open even when wearing gloves

For more information or to request a sample, e-mail infoeu@contecinc.com or phone + 33 (0) 2 97 43 76 98.

Copyright© 2019 Contec, Inc. All rights reserved. INT120 181019

Contec
Vannes Cedex
France
Tel: +33(0) 2 97 43 76 98

Contec Inc
Spartanburg, SC
United States
Tel: +1 (864) 503-8333

Contec Cleanroom Technology (Suzhou) Co. Ltd
Suzhou
China
Tel: +86-512-6274 4050



Product Specification

Contec HydroPure 500ml Capped Bottle

6% hydrogen peroxide in purified water
 Ready to use with a 12 month shelf life
 0.2 micron filtered
 Aseptic filled under UDAF grade A in a grade B cleanroom
 Double bagged in polyethylene linear tear bags
 Relevant certificates of analysis and sterility provided with every batch.
 Full supporting documentation information is available on request.

SATWIPE Canister 100 wipes

54gsm polyester/cellulose wipes 150 x 230mm
 12 month shelf life
 Packed in a grade C cleanroom
 Gamma irradiated at no less than 25kGy

Efficacy

Contec HydroPure has been showed to have good efficacy against spores against standard laboratory tests. Contec HydroPure SATWIPES have also been tested against the EN 16615 test for presaturated wipes against bacteria, fungi and bacterial spores.

Test	Test Type	Organism	Contact Time	Result	Test	Test Type	Organism	Contact Time	Result*
EN1276	Suspension Test	<i>E. coli</i>	5 min	Log 5.3	EN16615	Wiping Test	<i>S. aureus</i>	2 min	Log 5.33
EN1650	Suspension Test	<i>A. brasiliensis</i>	15 min	Log 4.0	EN16615	Wiping Test	<i>P. aeruginosa</i>	2 min	Log 5.34
EN13704	Suspension Test	<i>B. subtilis</i>	15 min	Log 3.0	EN16615	Wiping Test	<i>A. brasiliensis</i>	2 min	Log 4.29
EN13697	Surface Test	<i>B. subtilis</i>	15 min	Log 3.0	EN16615	Wiping Test	<i>B. subtilis</i>	2 min	Log 3.18

Product Application

Pour 500ml Contec HydroPure into the centre of the wipe roll. Close lid and invert canister multiple times until all fluid is absorbed into the wipes, this should take less than 5 mins. Once saturated with Contec HydroPure, the wipes are suitable for sporicidal disinfection. Wipes should be used within 7 days and then disposed of. There is space to write the saturation date on the front of the canister.

Part No.	Description	Size	Packaging
SSWC16HP	Contec Sterile Amplitude Theta SATWIPES System For use with Contec Sterile HydroPure	150 x 230mm	100 wipes per canister; 10 canisters per case

Use biocides safely. Always read the label and product information before use.