

SATWIPES® Sigma™ Wipes

Cellulose/polyester wipes presaturated with IPA/DI water


SATWIPES® Sigma™ consists of a perforated roll of saturated wipes which are dispensed through a center-pull canister. The wipes are made of 55% cellulose and 45% polyester creating highly sorbent and durable wipe. The isopropyl alcohol solvent blends provide safe cleaning solutions in various concentrations. Packaged with refill packs and a reusable pop-up dispenser, eliminating the need for solvent bottles and reduces solvent handling.



Features	Benefits
Hydraulically entangled fabric	<ul style="list-style-type: none"> • Creates a soft yet durable wipe
Clean wipe	<ul style="list-style-type: none"> • Offers consistent solvent application
Low levels of particles and extractable contaminants	<ul style="list-style-type: none"> • Reduces contamination possibility

Part No.	Description	Size	Packaging
SAT-C1-100	SATWIPES Sigma Wipes, Presaturated with 99% IPA	6" x 9" (152 x 230 mm)	100/roll; 12 rolls + canister/case
SWC10008	SATWIPES Sigma Wipes, Presaturated with 99% IPA	6" x 9" (152 x 230 mm)	100/roll; 10 rolls
SWC10190	SATWIPES Sigma Wipes, Presaturated with 99% IPA	6" x 9" (152 x 230 mm)	100/roll; 18 rolls
SAT-C1-7030	SATWIPES Sigma Wipes, Presaturated with 70% IPA/30% DI Water	6" x 9" (152 x 230 mm)	100/roll; 12 rolls + canister/case
SAT-C1-7030/18	SATWIPES Sigma Wipes, Presaturated with 70% IPA/30% DI Water	6" x 9" (152 x 230 mm)	100/roll; 18 rolls
SAT-C1-1090	SATWIPES Sigma Wipes, Presaturated with 70% IPA/30% DI Water	6" x 9" (152 x 230 mm)	100/roll; 12 rolls + canister/case

Product Information

Material	<ul style="list-style-type: none"> 55% cellulose/45% polyester
Construction	<ul style="list-style-type: none"> Hydroentangled nonwoven
Packaging materials	<ul style="list-style-type: none"> Case (CS), corrugated fiberboard (PAP) 



Recycle Symbols



Technical Data

	99% IPA	70/30% IPA	90/10% IPA	
Attribute (units)	Typical Value	Typical Value	Typical Value	Test Method
Basis weight, nominal; (g/m ²)	66	68	68	Contec Method
Non-volatile residue, NVR				IEST-RP-CC004.2, Sec. 6.1.2
In deionized water; (g/m ²)	0.056	0.056	0.056	
In isopropyl alcohol; (g/m ²)	0.015	0.015	0.015	
Particles, readily releasable				IEST-RP-CC004.2, Sec. 5.1
Particles ≥ 0.5µm; (x10 ⁶ /m ²)	20.0	17.1	17.1	
Fibers ≥ 100µm; (x 10 ³ /m ²)	23	14.8	14.8	IEST-RP-CC004.2, Sec. 5.2

Packaging	EA/RL	RL/PCH	PCH/CS	EA/CS
SAT-C1-100	100	1	12	1,200
SWC10008	100	1	10	1,000
SWC10190	100	1	18	1,800
SAT-C1-7030	100	1	12	1,200
SAT-C1-7030/18	100	1	18	1,800
SAT-C1-1090	100	1	12	1,200

EA = each; PCH = pouch; RL = roll; CS = case

Notes

- The data shown are typical values and should not be used as product specifications.
- Valid product comparisons may only be obtained through side-by-side testing in the same test facility, under similar conditions.
- Current and/or comparison data may be available. Please contact a Contec sales representative for details.
- All of Contec's packaging is compatible with hydrogen peroxide gassing applications.
- These wipes are free of lint and loose fibers, and meet the definition of lint-free/low linting wipes according to the United States Pharmacopoeia Chapter 797 (USP-NF General Chapter <797> Pharmaceutical Compounding -Sterile Preparations) and the Institute of Environmental Sciences and Technology Recommended Practice IEST-RP-CC004.4

Test Methods:

- CTM = Contec Test Method
- IEST-RP-CC004.3 = Evaluating Wiping Materials Used in Cleanroom and Other Controlled Environments, Institute of environmental Sciences and Technology, Rolling Meadows IL

PDSW050 | 091922

Copyright © 2022 Contec, Inc. All rights reserved.

Contact us:
info@contecinc.com
contecinc.com

Contec, Inc.
Spartanburg, SC 29303 USA
tel: +1-864-503-8333

