

Contec® Healthcare TB1-3300™ Disinfectant

Ethyl Alcohol/Quat Cleaner Disinfectant

Contec® Healthcare TB1-3300™ Disinfectant is a blend of quaternary ammonium and ethyl alcohol available in 32 oz. (0.9 L) and 1 gallon (3.8 L) containers. TB1-3300 Disinfectant is available sterile and nonsterile. It is a ready to use cleaner/disinfectant that generates low residue.

TB1-3300 Disinfectant has a 1-minute bactericidal, virucidal, tuberculocidal and fungicidal contact time. It is bleach-free, making it an excellent choice for operators in critical environments, sterile compounding facilities, retail pharmacies, and laboratories.

A flip-top cap is included for controlled application. Contec Healthcare TB1-3300 Disinfectant is ideal for cleaning and disinfecting environmental surfaces in pharmacy cleanrooms and laboratory areas including work stations, isolators, Primary Engineering Controls (PECs) (sterile TB1-3300), Secondary Engineering Controls (SECs), and other hard non-porous environmental surfaces.

Please contact your sales representative for further information.





Features	Benefits
One-step cleaner/disinfectant	<ul style="list-style-type: none"> • Compatible with common cleanroom surfaces • Great disinfectant/germicidal detergent for daily use
1-minute contact time	<ul style="list-style-type: none"> • Easy to maintain compliance and appropriate dwell time
Low residue	<ul style="list-style-type: none"> • Less residue left behind after use • Residues are easily removed with common solutions*
0.2 micron filtered	<ul style="list-style-type: none"> • Appropriate for use in critical environments • No filtering required before use
Flip-top cap	<ul style="list-style-type: none"> • Better application control of disinfectant
Linear tear packaging	<ul style="list-style-type: none"> • Easy to open even when wearing gloves • Facilitates transfer into more critical areas
Bleach-free	<ul style="list-style-type: none"> • Compatible with common controlled environment surfaces

*Recommended rinse agents are isopropyl alcohol, ethanol, or sterile water.

Product Specifications	
Odor	Alcohol
EPA-Registered	Yes
% Ethyl Alcohol	72.5%
% Didecyl Dimethyl Ammonium Chloride	0.33%
Detergency	Excellent
Density	0.86*
pH	3.7 at 20°C
Shelf Life Unopened (sterile and nonsterile) Opened (sterile and nonsterile)	2 years 6 months (chemical potency)
HMIS Hazard Rating:	
Health	2
Flammability	3
Reactivity	0
Personal Protection	X

*See SDS for more detailed information.

Part no.	Description	Size	Packaging
HCTB3300	Contec Healthcare TB1-3300 Disinfectant, single bagged*	32 oz. (0.9 L)	12/case
HCTB3300IR	Contec Healthcare TB1-3300 Disinfectant, sterile, double bagged* 	32 oz. (0.9 L)	12/case
HCTB3305	Contec Healthcare TB1-3300 Disinfectant, single bagged	1 gallon (3.8 L)	4/case
HCTB3305IR	Contec Healthcare TB1-3300 Disinfectant, sterile, double bagged 	1 gallon (3.8 L)	4/case

Contec® Healthcare TB1-3300™ Disinfectant

EPA-Registered Disinfectant Efficacy Claims

Bactericidal Claim – 1 minute

<i>Bordetella pertussis</i>
<i>Burkholderia cepacia</i>
<i>Campylobacter jejuni</i>
<i>Enterobacter aerogenes</i>
<i>Escherichia coli</i> O157:H7 (E. coli)
<i>Klebsiella pneumoniae</i>
<i>Legionella pneumophila</i>
<i>Listeria monocytogenes</i>
<i>Proteus vulgaris</i>
<i>Pseudomonas aeruginosa</i>
<i>Salmonella enterica</i>
<i>Shigella dysenteriae</i>
<i>Staphylococcus aureus</i>
<i>Streptococcus pyogenes</i> (Strep)
<i>Streptococcus pneumoniae</i>

Drug Resistant and Multi-Drug Resistant (MDR) Bacteria – 1 minute

<i>Acinetobacter baumannii</i> – Multi-Drug resistant (MDR)
Carbapenem resistant <i>Enterobacter cloacae</i> (CRE) (KPC) (MDR)
Carbapenem resistant <i>Enterobacter cloacae</i> (CRE) (New Delhi) (MDR)
<i>Enterococcus faecium</i> – Vancomycin resistant (VRE)
<i>Escherichia coli</i> – ESBL Producing (MDR)
<i>Klebsiella pneumoniae</i> – ESBL Producing (MDR)
<i>Staphylococcus aureus</i> – Methicillin resistant (MRSA)
<i>Staphylococcus aureus</i> – Vancomycin resistant (VRSA)
<i>Staphylococcus aureus</i> Panel – Vancomycin intermediate (VISA)

Fungicidal Claim – 1 minute

<i>Candida albicans</i>
<i>Trichophyton interdigitale</i>

Virucidal Claim – 1 minute

Adenovirus type 5
Avian Influenza A Virus
Cytomegalovirus
Enterovirus EV-D68
Hepatitis B Virus (HBV)
Hepatitis C Virus (HCV)
Herpes Simplex Virus type 1
Human Coronavirus
Human Immunodeficiency Virus (HIV-1)
Influenza Virus type A2
Influenza B
Norovirus
Poliovirus serotype 1 (PV-1)
Respiratory Syncytial Virus
Rhinovirus
Rotavirus Strain WA
Vaccinia Virus

Animal Pathogens – 1 minute

Avian Influenza A Virus
<i>Bordetella bronchiseptica</i>
Bovine Diarrhea Virus (BVDV)
Canine distemper Virus
Canine parvovirus
Duck Hepatitis B Virus
Feline calicivirus
Feline leukemia Virus

Tuberculocidal Claim – 1 minute

<i>Mycobacterium bovis</i>
