

PREempt[®] WIPES

One-Step Surface Cleaner and Disinfectant Wipes formulated

with **ACCELERATED[®]**
HYDROGEN PEROXIDE
Patented

WHY CHOOSE PREempt[™]?

- **CLEANER** Effective disinfectant as well as a cleaner resulting in added confidence that disinfection can occur.
- **FASTER** Rapid and realistic contact times on a broad spectrum of difficult to kill pathogens.
- **RESPONSIBLE** Designed to be easier on surfaces, staff and occupants. Users are not exposed to VOCs (volatile organic compounds) and PPE** is not required as per health and safety rating.
- **SUSTAINABLE** The active ingredient, Hydrogen Peroxide, breaks down into water and oxygen leaving no active residues which helps to reduce environmental impact.

DIRECTIONS FOR USE

Disinfection of non-porous hard surfaces such as the exterior of fume hoods, table tops and other environmental surfaces.

1. Heavily soiled surfaces require cleaning prior to disinfection.
2. Apply to surface with disposable wipe.
3. *Ensure surface remains wet for 3 minutes.*

Disinfection of Non-Critical, Pre-cleaned Instruments and Medical Devices:

1. Instruments must be thoroughly pre-cleaned to remove excess organic debris, rinsed, and dried.
2. Wipe the instrument with PREempt[™] Wipes and allow surface to remain wet for 3 minutes.
3. Wipe with a clean, damp cloth or paper towel and allow to air dry.

**PPE - Personal Protective Equipment



For use in laboratories, clean rooms and other critical environments that require cleaning and surface disinfection. This includes work stations, fume hoods, laboratory counter tops, equipment and other hard non-porous environmental surfaces.



AVAILABLE PRODUCTS

PREempt[™] One-Step Surface Cleaner and Disinfectant Wipes

DIN: 02436671

SHELF LIFE: 2 Years

CONCENTRATION OF ACTIVE:
0.5% Hydrogen Peroxide

Part Number: 11221

Size: 6" x 7" – 160ct wipes

Packaging: 12 canisters/case

PREempt™ One-Step Surface Cleaner Disinfectant Wipes

EFFECTIVE AGAINST A BROAD-SPECTRUM OF BACTERIA, VIRUSES AND FUNGI IN 3 MINUTES!

Virucidal: 3 Min.

- Proven effective against the Poliovirus Type 1, Sabin strain type 1 which allows for a Broad-Spectrum Virucidal claim against most enveloped and non-enveloped viruses
- HIV-1 Human Immunodeficiency Virus (HIV)
- Human Coronavirus 229E
- Feline Calicivirus, F9 Strain, as a surrogate for Norwalk and Norwalk-like viruses
- This product has demonstrated effectiveness against Poliovirus and is expected to inactivate all Influenza A viruses including 2009 (H1N1) pandemic Influenza A virus

Bactericidal: 3 Min.

- *Staphylococcus aureus*
- *Pseudomonas aeruginosa*
- *Salmonella enterica*
- Methicillin-resistant *Staphylococcus Aureus* (MRSA)
- Vancomycin-resistant *Enterococcus faecalis* (VRE)

Fungicidal: 3 Min.

- *Trichophyton mentagrophytes*

Tuberculocidal: 3 Min.

- *Mycobacterium terrae*

Broad-Spectrum Sanitizing: 30 Sec.

- *Staphylococcus aureus*
- *Pseudomonas aeruginosa*
- *Salmonella enterica*
- *Klebsiella pneumoniae*
- *Escherichia coli* 0157:H7
- *Campylobacter jejuni*

PREempt™ One-Step Surface Cleaner and Disinfectant Wipes are non-woven, meltblown polypropylene wipe material, ensuring even surface coverage and increasing the ability and efficiency of AHP® for fine particle removal.

*Please refer to reference sheet for a complete list of pathogenic organisms, additional features, and claims.

PREempt™ branded products are manufactured by Virox Technologies Inc.

Accelerated Hydrogen Peroxide® (AHP®) and PREempt™ are registered trademarks of Virox Technologies Inc.



Manufactured by



*Engineering Revolutionary
Disinfectants for the War
Against Microbes*

Contec, Inc.
Spartanburg, SC 29303 USA
Toll Free: 1-800-289-5762 • +1-864-503-8333
healthcare@contecinc.com
www.contechealthcare.com