



PROSAT Theta Wipes

Cellulose/polyester nonwoven wipes presaturated with 70% IPA and deionised water

PROSAT Theta wipes are manufactured from 54g/m² lightweight hydroentangled cellulose/polyester and saturated with 70% IPA and deionised water.

Cellulose/polyester wipes are a cost-effective cleanroom wipe, with low levels of particles and fibres. Highly sorbent with good wet strength, the wipes can be used for many general cleaning applications. Ideal for wiping articles prior to pass through, routine cleaning and wipe down of lab tools, instruments and other equipment.

The wipes are provided in convenient and easy-to-use peel and reseal pouches. The resealable pouches preserve solvent saturation and cleanliness. PROSAT Sterile Theta wipes are gamma irradiated and validated per AAMI Guidelines to a 10⁻⁶ SAL, so can be used in all grades of life science cleanroom.

When used as a disinfectant, the IPA wipes are efficacious against bacteria in 1 min and yeasts in 3 mins. PROSAT Theta wipes are authorised for sale in the EU and United Kingdom under the EU and GB Biocidal Products Regulation.



Features	Benefits
Hydroentangled cellulose/polyester fabric	<ul style="list-style-type: none"> Highly sorbent with good wet strength Low in particles and fibres Excellent general purpose wipe
No binders or additives	<ul style="list-style-type: none"> Leaves no residue on the cleanroom surface
Resealable pouch	<ul style="list-style-type: none"> Maintains saturation levels throughout use and makes wipe removal easy
Presaturated wipes	<ul style="list-style-type: none"> Reduces solvent usage and VOC emissions Ensures consistent saturation of each wipe independent of the operator
Validated sterile via gamma irradiation to a 10 ⁻⁶ Sterility Assurance Level	<ul style="list-style-type: none"> Suitable for use in Grade A/B cleanrooms

Part No.	Description	Size	Packaging
PSC20005-BPR	PROSAT Sterile Theta Wipes Presaturated with 70% IPA and 30% DI Water	230 x 280mm	40 wipes per pouch 40 pouches per case
PSC20006-BPR	PROSAT Theta Wipes Presaturated with 70% IPA and 30% DI Water	230 x 280mm	50 wipes per pouch 28 pouches per case

For more information or to request a sample, email infoeu@contecinc.com or phone +33 (0) 2 97 43 76 98.

Copyright © 2023 Contec, Inc. All rights reserved. INT190 130623



Product Information	
Material and construction	Hydroentangled 54% cellulose / 46% polyester
Saturant	70% IPA (USP Grade) with 30% DI water
Sterility	Gamma irradiated to a validated SAL 10 ⁻⁶
Shelf life	Sterile: 2 years from manufacturing date
Environment	ISO 5 - 8 Grade A/B for sterile, C/D for nonsterile

Technical Data		
Attribute (units)	Typical Value	Test Method
Basis weight, nominal; (g/m ²)	54	Contec Method
Non-volatile residue, NVR		IEST-RP-CC004.3, Sec. 7.1.2
In deionized water; (g/m ²)	0.037	
In isopropyl alcohol; (g/m ²)	ND	
Specific ions		IEST-RP-CC004.3, Sec. 7.2.2
Sodium; (ppm)	11.9	
Chloride; (ppm)	18.2	
Particles, readily releasable		
Particles ≥ 0.5µm; (x10 ⁶ /m ²)	18.6	IEST-RP-CC004.2, Sec. 5.1
Fibres ≥ 100µm; (x 10 ³ /m ²)	21.7	IEST-RP-CC004.2, Sec. 5.2

Notes

- a) The data shown are typical values and should not be used as product specifications.
- b) Valid product comparisons may only be obtained through side-by-side testing in the same test facility, under similar conditions.
- c) Current and/or comparison data may be available. Please contact a Contec sales representative for details.

Test Methods

CTM Contec Test Method
 IEST-RP-CC004.3 Evaluating Wiping Materials Used in Cleanroom and Other Controlled Environments, Institute of Environmental Sciences and Technology, Rolling Meadows IL



Recycling Key



Efficacy Information

Test	Description	Log Reduction	Time	Test	Description	Log Reduction	Time
EN16615	<i>E. hirae</i>	>5.03	1 min	EN16615	<i>P. aeruginosa</i>	>5.09	1 min
EN16615	<i>S. aureus</i>	>5.32	1 min	EN16615	<i>C. albicans</i>	>4.06	3 min

Packaging Information

Packaging Materials						VOC Content		
	Outer bags	Low density polyethylene (LDPE) 				kg/CS	kg/PCH	
	Case	Corrugated fibreboard (PAP) 						
Packaging Configuration		EA/PCH	PCH/OB1	OB1/OB2	OB2/CS	EA/CS		
PSC20005-BPR		50	1*	7	4	1400	8.45	0.30
PSC20006-BPR		50	4	7		1400	8.96	0.32

EA = Each, PCH = Pouch, OB = Outer Bag 1/2/3, CS = Case
 *Not suitable for recycling

Use biocides safely. Always read the label and product information before use.