



CONSTIX[®] Swab, SV-2

Sealed polyester swabs for surface cleaning validation

Description

Designed for validation of surface cleaning protocols in bio-pharmaceutical industries, the SV-2 package contains five SP-1 polyester swabs specially processed to minimize background TOC contribution from the swab.

The large surface of the SP-1 swab head and the long flexible handle provide exceptional surface coverage when swabbing. This ensures the entire surface of interest is swabbed effectively and maximizes retrieval.

Features and Benefits

- Low background TOC content for minimal interference with test results
- Strong, durable polyester swabs for optimum retrieval
- Sterilizable by most common methods
- Five swabs per pack allows complete usage with no waste

Applicability

The information presented here is applicable to the part numbers shown below, as well as to any product containing the same materials and produced under the same conditions, regardless of product size or packaging configuration. Please contact a Contec sales representative for more details.

SV-2

Technical Data

Swab head material	100% knitted polyester
Swab handle material	100% virgin polypropylene
Construction	Ultrasonically sealed

Attribute; (units)	Value **	Test Method
Head dimensions		
Width; (in./mm)	0.25/6.0	
Length; (in./mm)	0.8/20.3	
Handle dimensions		
Width; (in./mm)	0.12/3.0	
Thickness; (in./mm)	0.12/3.0	
Overall length; (in./mm)	6.0/152.4	
Total Organic Carbon (TOC); (ppm, per swab)	< 0.15	

** ND = None detected; levels are below detection limit of test equipment

Notes:

- 1) Data shown are typical values and should not be used as product specifications.
- 2) Valid product comparisons may only be obtained through side-by-side testing in the same test facility, under similar conditions.
- 3) Current and/or comparison data may be available. Please contact a Contec sales representative for details.

Revision date: 09/05/06



P.O. Box 530
 Spartanburg, SC 29304 USA
 Phone: +1-864-503-8333
 Fax: +1-864-503-8444
 Toll free US, Canada, PR: 1-800-289-5762

Customer Service
 United States: 1-866-855-4682
 Europe/Asia: +1-864-699-8257
 e-mail: wipers@contecinc.com
 www.contecinc.com