



TuffStuff™ Creped White Wipes

Critical Task Wipes

Description

TuffStuff Creped Wipes are white nonwoven creped wipes made from 55% cellulose and 45% polyester. They offer low linting, strength, good absorbency and versatility. Ideal for a wide range of critical cleaning tasks essential to industries such as biotech, pharmaceutical, and electronics. They are very effective in applications such as surface preparations for coatings, machine operations and general maintenance. Tuff Stuff wipes are available in quarterfolded and dispensing roll configurations.

Features and Benefits

- Clean, contains no binders
- Very sorbent
- Low linting
- Durable maintains strength when wet
- Solvent resistant
- Creped for softness and sorbency

Applicability

The information presented here is applicable to the part numbers shown below as well as to any product containing the same materials and produced under the same conditions, regardless of product size or packaging configuration. Please contact a Contec sales representative for more details.

TSWC1214, CIC1217Q

Technical Data

Attribute; (units)	Value **	Test Method
Basis weight; (g/m ²)	71	ASTMD3776
Thickness (mils)	19	ASTM D1777
Sheet Grab Tensile (lbs)		
Machine Direction (MD)	33	ASTM D5034
Cross Direction (XD)	20	ASTM D5034
Sorbency in water		
Intrinsic; (mL/g)	4.9	ASTM D6651-01
Extrinsic; (mL/m ²)	345	ASTM D6651-01
Sorptive rate, (seconds)	< 1	

** ND = None detected; levels are below detection limit of test equipment
Prelim 7/10/2007

Notes:

- 1) Data shown are typical values and should not be used as product specifications.
- 2) Valid product comparisons may only be obtained through side-by-side testing in the same test facility, under similar conditions.
- 3) Current and/or comparison data may be available. Please contact a Contec sales representative for details.

Revision date: 02/18/09
ETR 1490



P.O. Box 530
Spartanburg, SC 29304 USA
Phone: +1-864-503-8333
Fax: +1-864-503-8444
Toll free US, Canada, PR: 1-800-289-5762

Customer Service
United States: 1-866-855-4682
Europe/Asia: +1-864-699-8257
e-mail: wipers@contecinc.com
www.contecinc.com