

# Contec® Sterile 70% Isopropanol

## USP <797> Compliant Sterile Alcohol

Contec Sterile 70% Isopropanol contains 70% by volume USP Grade Isopropanol (IPA) and 30% USP Grade Purified Water. Filtered to 0.2µ and induction sealed to assure sterility, the bottles are double-bagged and gamma-irradiated. Each bottle has a lot code and expiration date for easy record keeping. Validated per the Association for the Advancement of Medical Instrumentation (AAMI) 11137:2006 Guidelines to a 10<sup>-6</sup> Sterility Assurance Level (SAL). *Note: Sterility validation is performed after trigger is installed.*

A convenient trigger sprayer with settings for stream or mist spraying is included. The stream setting is ideal for targeting specific hard to reach areas, while the mist setting is used to cover larger areas. Each leak-proof bottle is induction sealed, capped and shipped with an uninstalled trigger sprayer to eliminate the occurrence of leakage. Instructions for sterile installation of the trigger sprayer are available.

### Applications:

- Wipe down of isolators, cabinets, laminar flow hoods and surfaces
- Use on gloved hands in sterile suites
- Wipe down of products and equipment for pass through to controlled environments
- Ideally suited for use with sterile cleanroom wipes



<b>Appearance</b>	Clear, colorless liquid
<b>Odor</b>	Characteristic alcohol
<b>Water Quality</b>	Meets the requirements of USP Grade Deionized Water
<b>% IPA Contents (+/-3%)</b>	67.9% - 72.1%
<b>Specific Gravity (20°C)</b>	0.857 - 0.867

FEATURE	BENEFIT
Contains 70% USP-grade Isopropanol	<ul style="list-style-type: none"> <li>• Compliant with USP &lt;797&gt;</li> <li>• Evaporates and leaves extremely low residues</li> <li>• Eliminates mixing, filtration, processing and QC expense of in-house production of sterile Isopropanol water blends</li> </ul>
Filtered through 0.2 micron in-line filter	<ul style="list-style-type: none"> <li>• Eliminates the need for in-house filtration</li> <li>• Meets USP Grade Purified Water</li> </ul>
Gamma-irradiated to a sterility assurance level (SAL) of 10 <sup>-6</sup>	<ul style="list-style-type: none"> <li>• Compliant with USP &lt;797&gt;</li> <li>• Prefiltered, packaged and sterilized for easy introduction into critical environments</li> <li>• Certificate of processing and sterility with each case</li> </ul>
Induction sealed bottles	<ul style="list-style-type: none"> <li>• Reduces opportunity for leakage during sterilization and shipment</li> </ul>
Adjustable trigger spray	<ul style="list-style-type: none"> <li>• Select spray or jet option</li> </ul>
Plastic bottle	<ul style="list-style-type: none"> <li>• Recyclable, environmentally friendly, no costly disposal</li> </ul>
Double bagged	<ul style="list-style-type: none"> <li>• To preserve product integrity</li> </ul>
Lot numbered bottles	<ul style="list-style-type: none"> <li>• Assures complete product traceability</li> </ul>

NOT FOR RETAIL SALE. MATERIAL DATA SHEET AVAILABLE FOR PROFESSIONAL USER UPON REQUEST.

Part No.	Description	Quantity
SB1670301R	Sterile 70% Isopropanol, 16 oz. (0.47L) w/trigger sprayer	12 bottles/case
SB3270301R	Sterile 70% Isopropanol, 32 oz. (0.95L) w/trigger sprayer	12 bottles/case
SB12870301R	Sterile 70% Isopropanol, 1 gallon (3.8L) bottle	4 bottles/case



P.O. Box 530 • Spartanburg, SC 29304 USA  
 Tel: +1-864-503-8333  
 Toll Free (US, Canada, PR): 1-800-289-5762  
 Europe: +33 297437690  
 www.contecinc.com • disinfect@contecinc.com