

Critical Environments Mop Guide

Brand Name	Style	Composition	Options Available		Operational Characteristics				Areas of Use						
			Presaturated	Sterile	Autoclavable	Low Particle Generation	Absorbency	Abrasion Resistance	Residue Removal	Disinfectant Application	Wet Cleaning	Dry Cleaning	Floors	Walls & Ceilings	RABS & Isolators
Edgeless®	Tubular loop	100% knitted polyester	No	Yes	Yes	Excellent	Good	Excellent	Excellent	Excellent	Excellent	Poor	Yes	No	No
VertiKlean®	Sponge	Ester foam laminated to polyester fabric	Yes	Yes	Yes	Good	Good	Good	Good (textured)	Excellent	Good	Fair	No	Yes	No
VertiKlean® MAX™	Sponge	P U foam laminated to polyester fabric	Yes	Yes	No	Good	Excellent	Good	Fair	Excellent	Good	Fair	Yes	Yes	Yes
EasyReach™	Flat mop	2-ply 100% knitted polyester	Yes	Yes	No	Excellent	Good	Good	Fair	Excellent	Good	Fair	No	No	Yes
QuickTask® 0200	Flat mop	100% knitted polyester	No	Yes	No	Excellent	Good	Good	Fair	Excellent	Good	Fair	Yes	Yes	No
QuickTask® 0300	Flat mop	Microfiber	No	No	Yes	Fair	Excellent	Good	Good	Good	Good	Good	Yes	Yes	No
QuickTask® 0400	Flat mop	Microfiber	No	No	Yes	Good	Excellent	Good	Good	Good	Good	Good	Yes	Yes	No
QuickTask® 0500	Flat mop	Polyester/cellulose blend	No	Yes	No	Fair	Good	Fair	Excellent	Good	Good	Good	Yes	Yes	No
EasyCurve™	Sponge	P U foam laminated to polyester fabric	No	Yes	No	Good	Excellent	Good	Fair	Excellent	Good	Fair	Yes	Yes	No
EasySat™	Sponge	Ester foam laminated to polyester fabric	No	Yes	Yes	Good	N/A	Good	Fair	Excellent	Good	N/A	Yes	No	No
Tax-Fre®	Dry Wipe	Low-tack treated nonwoven polyester	No	No	No	Good	N/A	Good	N/A	N/A	N/A	Excellent	Yes	Yes	No