

# SterileSorb™ Sterile Nonwoven Wipe

Contec® introduces SterileSorb™ sterile hydroentangled polyester/cellulose wipe

Contec's SterileSorb™ sterile nonwoven wipes are ideal for general purpose cleaning in sterile environments. Packaged in our sterile peelable outer bag and linear tear inner bag specifically designed for sterile products. Embossed lot numbers on both inner and outer bags ensure lot traceability. Sterile validation reports are available upon request.

## Features

- Comprised of strong, durable hydroentangled fabric
- Sterilized using gamma irradiation
- Extremely absorbent
- Low levels of particles and extractables
- Excellent for wiping and absorbing spills in cleanroom environments
- Available in 9" x 9" and 12" x 12" sizes



*SterileSorb™ nonwoven wipes are packaged in our sterile peelable outer bag and linear tear inner bags.*

Part #	Description	Piece/IB	IB/OB	OB/Case
C2-99IR/25	9" x 9" (22.9 cm x 22.9cm) wipe	25	6	12
C2-1212IR	12" x 12" (30.5 cm x 30.5cm) wipe	25	6	7

Note: IB= Inner Bags, OB= OuterBags



 **CONTEC®**  
An ISO 9001:2008 registered company

P.O. Box 530 • Spartanburg, SC 29304 USA  
+1-864-503-8333 • fax +1-864-503-8444  
Toll free US, Canada, PR 1-800-289-5762  
www.contecinc.com • e-mail: wipers@contecinc.com

Contec®, the Contec® logo and logotype, and SterileSorb™ are registered trademarks of Contec, Inc. Contec reserves the right to make changes to this sheet. ©2010 Contec, Inc.

WIP002  
091907