



Contec® Irradiated 70% IPA LE

Low endotoxin irradiated alcohol with Water for Injection

Contec Low Endotoxin Irradiated Alcohol contains 70% by volume USP Grade Isopropanol (IPA) and 30% Water for Injection (WFI). The alcohol blend is 0.2 micron filtered, filled, and bagged in an ISO Class 7 cleanroom. This clean manufacture coupled with water for injection means the alcohol blend is guaranteed and certified to have an endotoxin level < 0.25 EU/ml, making it suitable for product contact surfaces. The bottles are provided with the triggers installed and gamma irradiated at no less than 25kGy. Each lot is sterility tested and certified.

Available as a 17 oz and 34 oz trigger spray bottle, the trigger sprays are fitted with a protected system, which ensure sterility throughout use. Each bottle has a lot code and expiration date for easy record keeping. The product has a 3 year unopened shelf-life from date of manufacture and a 6 month in-use shelf-life.



Trigger installed “bag in bottle” system protects contents during use.

FEATURES	BENEFITS
70% Isopropanol and 30% Water for Injection	<ul style="list-style-type: none"> • No residue for product contact areas
Guaranteed endotoxin level of less than 0.25 EU/ml	<ul style="list-style-type: none"> • Ideal for use in product contact areas
Double bagged in linear tear packaging	<ul style="list-style-type: none"> • Easy to open even when wearing gloves • Facilitates transfer into isolators and RABS
Trigger installed “bag in bottle” protected system	<ul style="list-style-type: none"> • Eliminates the risk of the alcohol being contaminated when installing the trigger • Prevents bottle contents from being contaminated during use
Gamma irradiated at no less than 25kGy and lot tested for sterility	<ul style="list-style-type: none"> • Suitable for use in the most critical environments
Lot numbered bottles	<ul style="list-style-type: none"> • Assures complete product traceability
Trigger spray can be set to jet or spray	<ul style="list-style-type: none"> • Large droplet size reduces the risk of inhalation and provides good surface coverage

Part No.	Description	Quantity
SBT167030LE	Contec Irradiated 70% IPA/30% WFI LE, 17 oz. (0.5L), trigger installed, sterile	8 bottles/case
SBT347030LE	Contec Irradiated 70% IPA/30% WFI LE, 34 oz. (1L), trigger installed, sterile	6 bottles/case

Visit our website for more information: www.contecinc.com.

Copyright © 2018 Contec, Inc. All rights reserved. DIS034 101518

Contec, Inc
Spartanburg, SC
United States
Tel: +1 (864) 503-8333

Contec
Vannes Cedex
France
Tel: +33(0) 2 97 43 76 98

Contec Cleanroom Technology (Suzhou) Co. Ltd
Suzhou
China
Tel: +86-512-6274 4050

