



Contec® HydroKlean

Sterile fill low endotoxin cleaner with hydrogen peroxide

Contec HydroKlean is a blend of 6% hydrogen peroxide and water for injection in a 17 oz and 34 oz (1L) bag in bottle with installed trigger sprayer and 170 oz (5L) capped container. It is ready to use and generates no Volatile Organic Compounds (VOC's).

Hydrogen peroxide breaks down to water and oxygen in use, so it leaves little to no residue on a surface. Contec HydroKlean has a certified endotoxin level of <0.25 EU/ml making it ideal for use in product contact areas.

The product is 0.2µm filtered and sterile filled into pre-irradiated containers under Grade A unidirectional airflow and bagged in a Grade B/ISO Class 5 cleanroom. The 17 oz (0.5L) and 34 oz (1L) "bag in bottle" system protects contents during use. Certificates of analysis are available for each batch.

6% hydrogen peroxide is not classed as corrosive and can be used safely in all areas of a cleanroom. It is nonflammable so it's ideal for use across large areas of a cleanroom or where sparks may be an issue. As Contec HydroKlean is very low residue, it eliminates the need for a residue removal stage with alcohol or water.

Contec HydroKlean is double or triple bagged allowing for ease of entry into controlled environments. The linear tear bags are easy to open even when wearing gloves.



HydroKlean "bag in bottle" system protects contents during use.

FEATURES	BENEFITS
Hydrogen peroxide breaks down to water and oxygen	<ul style="list-style-type: none"> • Very low residue and suitable for product contact areas • Nonflammable and generates no Volatile Organic Compound's (VOC's)
Guaranteed endotoxin level of less than 0.25 EU/ml	<ul style="list-style-type: none"> • Ideal for use in product contact areas
Sterile filled into pre-irradiated containers in an ISO Class 5 cleanroom	<ul style="list-style-type: none"> • Suitable for use in all areas of the cleanroom including ISO 5 cabinets and RABS
Double or triple bagged in linear tear packaging	<ul style="list-style-type: none"> • Easy to open even when wearing gloves • Facilitates transfer into isolators and RABS
Trigger spray in "bag in bottle" protected system	<ul style="list-style-type: none"> • Prevents bottle contents from being contaminated during use
Not classified as corrosive	<ul style="list-style-type: none"> • Can be used safely in all areas of the cleanroom

Part No.	Description	Quantity
SBT17HKLE	Contec HydroKlean, sterile fill, 17 oz (0.5L), trigger spray, triple bagged	8 bottles/case
SBT34HK6IR	Contec HydroKlean, sterile fill, 34 oz (1L), trigger spray, triple bagged	6 bottles/case
SBC170HKLE	Contec HydroKlean, sterile fill, 170 oz (5L), capped container, double bagged	2 containers/case

Visit our website for more information: www.contecinc.com.

Copyright © 2018 Contec, Inc. All rights reserved. DIS033 013018

Contec, Inc
Spartanburg, SC
United States
Tel: +1 (864) 503-8333

Contec
Vannes Cedex
France
Tel: +33(0) 2 97 43 76 98

Contec Cleanroom Technology (Suzhou) Co. Ltd
Suzhou
China
Tel: +86-512-6274 4050

