



PROSAT® Sterile™ Low Endotoxin Wipes

Polynit Heatseal wipes with 70% USP grade IPA,
30% Water for Injection

PROSAT Sterile Low Endotoxin Wipes are made 100% knitted polyester Polynit Heatseal wipes presaturated with a blend of 70% IPA and Water for Injection (WFI). Each batch of wipes has a guaranteed endotoxin level of less than <1EU/device. Suitable for use in aseptic suites.

The sealed edge wipe produces a very low level of particles and fibers and is ideal for use near product contact surfaces. Hanging slot in pouch makes it easy to use in isolators.

Each pouch contains a small number of wipes eliminating any waste at the end of a session. Double or triple bagged for easy entry into the sterile suites. Validated sterile, these Low Endotoxin wipes are suitable for use in the most critical pharmaceutical cleanrooms.



FEATURES	BENEFITS
Cleanroom laundered 100% knitted polyester wipes	<ul style="list-style-type: none"> • Very low levels of particles and fibers • Good abrasion and chemical resistance
Laser sealed edges	<ul style="list-style-type: none"> • Reduces fiber generation from the edge of the wipe
Presaturated with 70% IPA/WFI	<ul style="list-style-type: none"> • Controlled saturation levels greatly reduce solvent usage, VOC's and waste
Guaranteed endotoxin levels of <1EU/device	<ul style="list-style-type: none"> • Suitable for use near product contact surfaces
Validated sterile via gamma irradiation to a 10 ⁻⁶ Sterility Assurance Level	<ul style="list-style-type: none"> • Suitable for use in Grade A and B cleanrooms
Small number of wipes per resealable pouch	<ul style="list-style-type: none"> • Low number of wipes eliminates waste • Pouch can be easily resealed throughout use, maintaining saturation level
Double or triple bagged in linear tear packaging	<ul style="list-style-type: none"> • Each bag is easy to open even when wearing gloves • Facilitates cleaning products entering isolators and RABS

Part No.	Description	Quantity
PSWE0001	PROSAT Sterile Polynit Heatseal with 70% IPA/30%WFI, 9" x 9" (23 x 23cm), triple bagged	10 wipes/pouch; 55 pouches/case
PSWE0002	PROSAT Sterile Polynit Heatseal with 70% IPA/30%WFI, 12" x 12" (30.5 x 30.5cm), double bagged	30 wipes/pouch; 15 pouches/case
PSWE0003	PROSAT Sterile Polynit Heatseal with 70% IPA/30%WFI, 12" x 12" (30.5 x 30.5cm), triple bagged	10 wipes/pouch; 30 pouches/case

Visit our website for more information: www.contecinc.com.

Copyright © 2018 Contec, Inc. All rights reserved. PRE005 101218

Contec, Inc
Spartanburg, SC
United States
Tel: +1 (864) 503-8333

Contec
Vannes Cedex
France
Tel: +33(0) 2 97 43 76 98

Contec Cleanroom Technology (Suzhou) Co. Ltd
Suzhou
China
Tel: +86-512-6274 4050

