



Contec® Sterile IPA

Contec® Sterile Isopropanol is a ready-to-use IPA solution pre-filtered to 0.2 micron and gamma irradiated. Containing USP Grade IPA and USP Grade Purified Water where applicable, it is designed for use in Pharmaceutical, Medical Device, and other Life Science industries. Each bottle is double-bagged to ensure sterility and printed with lot number and date of expiration for easy record keeping. Validated sterile to a 10⁻⁶ SAL per ANSI/AAMI/ISO 11137 guidelines. *Note: Sterility validation is performed after trigger sprayer is installed.*

A convenient trigger sprayer with settings for stream or mist spraying is included with 16 and 32 oz. bottles. The stream setting is ideal for targeting specific hard to reach areas, while the mist setting is used to cover larger areas. Each leakproof bottle is induction sealed, capped and shipped with an uninstalled trigger sprayer to eliminate the occurrence of leakage. Instructions for sterile installation of the trigger sprayer are available.



Product Application

- Wipe down of isolators, cabinets, laminar flow hoods and surfaces
- Use on gloved hands in sterile suites
- Wipe down products and equipment for pass through to controlled environments
- Ideally suited for use with sterile cleanroom wipe

Features	Benefits
Contains 70% USP grade isopropanol	<ul style="list-style-type: none"> • Compliant with USP <797> and USP <800> • Evaporates and leaves extremely low residues • Eliminates mixing, filtration, processing and QC expense of in-house production of sterile isopropanol water blends
Filtered through 0.2 micron in-line filter	<ul style="list-style-type: none"> • Eliminates the need for in-house filtration
Gamma irradiated to a sterility assurance level (SAL) of 10 ⁻⁶	<ul style="list-style-type: none"> • Prefiltered, packaged and sterilized for easy introduction into critical environments • Certificate of processing and sterility with each case
Plastic bottle	<ul style="list-style-type: none"> • Recyclable, environmentally friendly, no costly disposal

Part No.	Description	Packaging
SB167030IR	Sterile 70% Isopropanol, 16 oz. (0.47L), w/trigger sprayer	12/case
SB327030IR	Sterile 70% Isopropanol, 32 oz. (0.95L), w/trigger sprayer	12/case
SB1287030IR	Sterile 70% Isopropanol, 1 gallon (3.8L) bottle	4/case
SB128100IR	Contec Sterile Isopropanol, 128 oz. (3.8L)	4/case

Product Information

Appearance	Clear, colorless liquid
Odor	Characteristic alcohol
Water Quality	Meets the requirements of USP grade deionized water
% IPA Contents	70% 68.0% - 72.0%
	100% 99.0% - 99.99%
Specific Gravity (20°C)	70% 0.872 - 0.883
	100% 0.7861 - 0.7870
ISO Class cGMP Grade	ISO 3-8 Grade A/B
Shelf Life	Unopened 5 years
	Opened 6 months


Contec® Disinfectant and Solution Rotation

DAILY

Use everyday on all surfaces.

Disinfectants

- CyQuanol™
- Sterile CyQuanol™
- PREempt™ RTU
- PREempt™ Concentrate




ROUTINELY

Use based on room environment monitoring.

Sporicide

- PeridoxRTU®
- Sterile PeridoxRTU®
- Peridox® Concentrate




PERIODICALLY

Use to remove built up residue or when needed.

Detergent

- NeutraKlean
- Sterile NeutraKlean



RINSE

Use after cleaning and disinfecting.

Rinse

- IPA
- HydroKlean
- WFI
- Water

