


# Amplitude™ Epsilon™ Wipes

Amplitude™ Epsilon™ is made of a special grade blue cellulose and polyester hydroentangled to form a wipe with improved cleanliness and excellent sorbency. The cellulose uses a colorfast dye, to minimize color bleeding. Blue color allows for differentiation in critical areas. Good for wet areas requiring spill pickup and cleanup.



Features	Benefits
Colorfast blue dye	<ul style="list-style-type: none"> <li>• Allows for differentiation between wipes used in different areas</li> </ul>
Special grade blue hydroentangled cellulose and polyester	<ul style="list-style-type: none"> <li>• Colorfast dyeing coloring minimizes the risk of color bleeding</li> <li>• Hydroentangled to form a wipe with improved</li> </ul>
Excellent sorbency	<ul style="list-style-type: none"> <li>• Good for wet areas requiring spill pickup and cleanup</li> </ul>

Part No.	Description	Packaging
AMEP0001	Amplitude Epsilon Wipes, 9" x 9" (23x23cm), Flat stacked	300/bag; 12 bags/case
AMEP0002	Amplitude Epsilon Wipes, 12" x 12" (30.5x30.5cm), Flat stacked	150/bag; 20 bags/case

