

Contec® 70% Ethanol

70% CHP grade Ethanol and 30% purified water

Contec® 70% Ethanol contains 70% CHP grade Ethanol and 30% purified water. Filtered to 0.2µ and induction sealed to assure sterility, the bottles are double, easy tear bagged. Each bottle has a lot code and expiration date for easy record keeping.

A convenient trigger sprayer with settings for stream or mist spraying is included. The stream setting is ideal for targeting specific hard to reach areas, while the mist setting is used to cover larger areas. Each leak-proof bottle is induction sealed, capped and shipped with an uninstalled trigger sprayer to eliminate the occurrence of leakage. Instructions for sterile installation of the trigger sprayer are available. The shelf life is 5 years.



Product Application

- Wipe down of isolators, cabinets, laminar flow hoods and surfaces
- Wipe down products and equipment for pass through to controlled environments
- Ideally suited for use with sterile cleanroom wipes

Features	Benefits
Contains 70% CHP grade Ethanol and 30% purified water	<ul style="list-style-type: none"> • Evaporates and leaves extremely low residues
Filtered through 0.2 micron in-line filter in class 100 laminar hood	<ul style="list-style-type: none"> • Eliminates mixing, filtration, processing and QC expense of In-house production of ethanol water blends • Ensures the product is free from contamination and particulates
Induction sealed bottles	<ul style="list-style-type: none"> • Reduces opportunity for leakage during sterilization and shipment
Adjustable trigger spray	<ul style="list-style-type: none"> • Select spray or jet option
Plastic bottle	<ul style="list-style-type: none"> • Environmentally friendly, no costly disposal
Double easy tear bagged	<ul style="list-style-type: none"> • To preserve product integrity, even gloves are easy to tear
Lot numbered bottles	<ul style="list-style-type: none"> • Assures complete product traceability

Part No.	Description	Size	Packaging
SB1670EP	70% Ethanol, with trigger sprayer	16 oz. (0.47L)	12/case
SB3270EP	70% Ethanol, with trigger sprayer	32 oz. (0.95L)	12/case
SB16970EP	70% Ethanol, bottle	169 oz. (5L)	4/case

Product Information

Appearance	Clear, colorless liquid
Odor	Characteristic alcohol
Water Quality	Meets the requirements of CHP grade purified water
% IPA Contents	68% - 72%
Specific Gravity (20°C)	0.877 - 0.886
Packaging materials	Outer bags (OB1, OB2), low density polyethylene (LDPE)  Case (CS), corrugated fiberboard (PAP) 
Environment	ISO 5-8 Grade C/D
Shelf life	Unopened 5 years Opened 28 days

Recycle Symbols

