

Peroxigard®

Surface cleaner and disinfectant

Powerful. Responsible. Compatible. Faster.

Peroxigard® is an accelerated hydrogen peroxide (AHP). AHP® is a patented synergistic blend of hydrogen peroxide at low concentrations combined with surfactants, wetting agents, and chelating agents. This synergy dramatically increases hydrogen peroxide’s germicidal efficacy and cleaning performance.

All ingredients in AHP are listed on the EPA inerts lists and the FDA Generally Regarded as Safe List (GRAS). The benefits of AHP have been evaluated and supported by over 55 peer reviewed studies. Peroxigard has become the leading cleaner/disinfectant in the lab animal sciences space. It is compatible with most surfaces, and it is faster than other disinfectants in the market.



Features	Benefits
AHP technology	<ul style="list-style-type: none"> • Fast acting disinfectant activity • Can be used as a one-step disinfectant cleaner which saves time
Lowest EPA toxicity rating	<ul style="list-style-type: none"> • Odor-free for good operator acceptability • Contains no VOCs • Does not require personal protection equipment (PPE)
Contains hydrogen peroxide	<ul style="list-style-type: none"> • Breaks down into water and oxygen so low residue • Nonhazardous so allows cost effective disposal
Compatible with a wide range of materials	<ul style="list-style-type: none"> • Reduces corrosion and wear on surfaces and equipment

Part No.	Description	Size	Packaging
29101-CONT	Peroxigard RTU Disinfectant Includes flip-top cap	32 oz (0.9 L)	12/case
29105-CONT	Peroxigard RTU Disinfectant	1 gallon (3.8 L)	4/case
29109-CONT	Peroxigard RTU Disinfectant	5 gallon (19 L)	1/case
29110-CONT	Peroxigard RTU Disinfectant	30 gallon (114 L)	1/case
29111-CONT	Peroxigard RTU Disinfectant	55 gallon (208 L)	1/case
VHH073	Peroxigard RTU Secondary use label	-	10/pack
29305-CONT	Peroxigard Concentrate Surface Cleaner and Disinfectant	1 gallon (3.8 L)	4/case
29309-CONT	Peroxigard Concentrate Surface Cleaner and Disinfectant with spout	5 gallon (19 L)	1/case
29311-CONT	Peroxigard Concentrate Surface Cleaner and Disinfectant	55 gallon (208 L)	1/case
VHH074	Peroxigard Concentrate Secondary container label	-	10/pack
29221-CONT	Peroxigard RTU Disinfectant Wipes	6" x 7" (152 x 178 mm)	12/case
AHP425	Precise® Test Strips for Peroxigard Concentrate	-	50/bottle