
Amplitude[™] Sigma[™] Wipes

Cellulose/polyester hydroentangled nonwoven wipes

The Amplitude[™] Sigma[™] wipe is made from Sontara[®] nonwoven fabric comprised of a hydroentangled blend of 55% cellulose and 45% polyester staple fibers.

The Sigma wipe is strong and durable, even when wet, and is low in extractable residues and particles. It has outstanding sorbency and is an excellent all-purpose wipe.

Available validated sterile as Contec[®] SterileSorb[™].



Features	Benefits
Hydroentangled cellulose/polyester	Provides low levels of particles and extractables
Variety of fibers and blends	Excellent for wiping and absorbing spillsAllows the wipe to be strong and durable both wet and dry

Part No.	Description	Size	Packaging
AMSI0001	Amplitude Sigma Wipe	9" x 9"	300/bag;
	Flat stacked	(230 x 230 mm)	12 bags/case
AMSI0002	Amplitude Sigma Wipe	12" x 12"	150/bag;
	Flat stacked	(305 x 305 mm)	18 bags/case
AMSI0003	Amplitude Sigma Wipe	18" x 18"	75/bag;
	Flat stacked	(460 x 460 mm)	10 bags/case

PDSW004 | 080223 | ETR278, 4366

Copyright © 2023 Contec, Inc. All rights reserved.


Contact us: info@contecinc.com cleanroom.contecinc.com Contec, Inc. Spartanburg, SC 29303 USA tel: +1-864-503-8333



Product Information		
Material	Cellulose/polyester	
Construction	Hydroentangled nonwoven	
Packaging Materials	 Outer bags (OB1, OB2), low density polyethylene (LDPE) Case (CS), corrugated fiberboard (PAP) 	
Environment	ISO 5-8 Grade C/D	

Technical Data

Attribute (units)	Typical Value	Test Method		
Basis weight, nominal; (g/m²)	68	Contec Method		
Sorbent capacity; (mL/m²)	314	IEST-RP-CC004.3, Sec. 8.1		
Sorptive rate; (seconds)	<1	IEST-RP-CC004.3, Sec. 8.2		
Non-volatile residue, NVR		IEST-RP-CC004.3, Sec. 7.1.2		
In deionized water; (g/m²)	0.009			
In isopropyl alcohol; (g/m²)	0.009			
Specific ions		IEST-RP-CC004.3, Sec. 7.2.2		
Sodium; (ppm)	13.1			
Chloride; (ppm)	9.5			
Particles, readily releasable		IEST-RP-CC004.2, Sec. 5.1		
Particles $\geq 0.5 \mu m$; (x10 ⁶ /m ²)	32.7			
Fibers <u>></u> 100µm; (x 10 ³ /m ²)	22.0			







EA = each; OB = outer bag; CS = case

Notes

a) The data shown are typical values and should not be used as product specifications.

b) Valid product comparisons may only be obtained through side-by-side testing in the same test facility, under similar conditions.

c) Current and/or comparison data may be available. Please contact a Contec sales representative for details.

d) All of Contec's packaging is compatible with hydrogen peroxide gassing applications.

e) These wipes are free of lint and loose fibers, and meet the definition of lint-free/low linting wipes according to the United States Pharmacopoeia Chapter 797 (USP-NF General Chapter <797> Pharmaceutical Compounding -Sterile Preparations) and the Institute of Environmental Sciences and Technology Recommended Practice IEST-RP-CC004.4

Test Methods:

1. CTM = Contec Test Method

2. IEST-RP-CC004.3 = Evaluating Wiping Materials Used in Cleanroom and Other Controlled Environments, Institute of environmental Sciences and Technology, Rolling Meadows IL.

PDSW004 | 080223 | ETR278, 4366

Copyright © 2023 Contec, Inc. All rights reserved.

Contact us: info@contecinc.com cleanroom.contecinc.com Contec, Inc. Spartanburg, SC 29303 USA tel: +1-864-503-8333

