

Contec Filtered Alcohol

70% IPA or Denatured Ethanol blended with purified water

Contec alcohols are a blend of 70% v/v Isopropanol or denatured ethanol with 30% purified water. The alcohol blend is 0.2 micron filtered, filled and bagged in a Grade C cleanroom.

Contec’s ethanol is only denatured with IPA creating a very low residue product. Supplied as either trigger sprays fitted with a protected system, which ensures sterility throughout use, or 5L capped containers for when larger volumes are required.

Contec’s EU BPR product dossier for all products containing 70% IPA has received Union Authorisation from the ECHA Biocidal Products Committee. All Contec’s IPA products are authorised for sale in all EU countries and the UK. The Authorisation Number for the product family is EU-0020460-0000.



Features	Benefit
GMP manufactured under Grade A air flow in a Grade C cleanroom	<ul style="list-style-type: none"> Ensures the alcohol, container and packaging are free from contamination and particulates
No denaturants other than IPA	<ul style="list-style-type: none"> Completely residue free
Trigger sprayer and “bag in bottle ” protected system	<ul style="list-style-type: none"> Prevents bottle contents from becoming contaminated during use No need to discard unused product during shelf life Bottles can be completely emptied eliminating wastage
Trigger sprayer can be set to jet or spray	<ul style="list-style-type: none"> Large droplet size reduces the risk of inhalation and provides good surface coverage
Double bagged packed in linear tear packaging	<ul style="list-style-type: none"> Each bag is easy to open even when wearing gloves Facilitates transfer disinfection into cleanroom

For more information or to request a sample, email infoeu@contecinc.com or phone +33 (0) 2 97 43 76 98.

Copyright© 2022 Contec, Inc. All rights reserved. INT119 281222

