

Contec® CritiGear™ HD Chemo Gown

USP <800> compliant gown for use when compounding hazardous drugs*

Contec® CritiGear™ HD Chemo Gowns with BloxTech™ fabric technology are designed to provide protection against exposure to a wide range of hazardous drug compounds*. BloxTech's ASTM D6978-05 tested fabric provides wearer comfort as well as protection. The disposable gown's design features generous sleeve length with knit cuffs and thumb loops to ensure secure cuff area and comfort when performing compounding activities.

Contec CritiGear HD Chemo Gowns are constructed with USP <800> required sealed seams that also feature highly visible blue tape to ensure the correct gown is used for the correct application. The gown's wrap around design, with hook and latch closures on the back, helps to provide complete coverage during compounding activities. The extra-long ties, coupled with multiple size combinations, ensure a comfortable, secure fit. The extra-long ties also allow for the gown to be secured by tying the knot to the wearer's side, out of the way of compounding activities.



Contec CritiGear HD Chemo Gowns feature BloxTech fabric that is proven to protect the wearer from a wide range of hazardous drugs*.

Features	Benefits
BloxTech fabric technology	<ul style="list-style-type: none"> • Low-lint disposable fabric is compliant with the 2022 revision to USP Chapter <797>
Tested per ASTM D6978-05 standards	<ul style="list-style-type: none"> • *Specifically tested against <i>Carmustine, Cisplatin, Cyclophosphamide, Doxorubicin, Etoposide, Fluorouracil, Methotrexate, Paclitaxel, and Thiotepa</i>
USP <800> compliant	<ul style="list-style-type: none"> • Constructed according to USP <800> requirements for hazardous drug compounding
Complete protection	<ul style="list-style-type: none"> • Full closure back, hook and loop neck closure, and taped seams provide total coverage against permeation
Thumb loop	<ul style="list-style-type: none"> • Knitted cuffs with thumb loops eliminates movement of sleeves
Easy removal design	<ul style="list-style-type: none"> • Hook and loop closure allows user to simply untie side tie and pull gown down and away from the body with one hand
Combination sizing	<ul style="list-style-type: none"> • Combination sizing for ease of donning over other protective frocks such as Contec CritiGear Cleanroom Frocks

Part no.	Description	Size	Packaging
HCGA0800-S/M	Contec CritiGear HD Chemo Gown	Small/Medium	10/bag, 3 bags/case
HCGA0800-L/XL	Contec CritiGear HD Chemo Gown	Large/X-Large	10/bag, 3 bags/case
HCGA0800-2XL/3XL	Contec CritiGear HD Chemo Gown	2X-Large/3X-Large	10/bag, 3 bags/case
HCGA0800-4XL/5XL	Contec CritiGear HD Chemo Gown	4X-Large/5X-Large	10/bag, 3 bags/case
HCGA0800-6XL/7XL	Contec CritiGear HD Chemo Gown	6X-Large/7X-Large	10/bag, 3 bags/case