

EasyCurve™ Mop

Lightweight mop for ceilings, walls and floors

Contec's EasyCurve™ Mop was designed specifically for critical and controlled environments, now available with a ribbed microfiber face giving users smooth and easy mopping for better particle pick up and contamination removal. The combination of the curved mop face and the pivoting connector results in effortless surface contact and outstanding maneuverability, making the EasyCurve an ideal cleaning tool for ceilings, walls and floors. Available in non-sterile and validated sterile versions.

All components of the EasyCurve Mop are autoclavable. For convenience, EasyCurve mop heads are available Validated Sterile.

Applications

Designed for use in pharmaceutical and medical device manufacturing plants, and other similar critical or controlled environments, the EasyCurve is an outstanding all-surface mop for use on ceilings, walls and floors.



With its absorbent and drip-free polyurethane foam mop head, the EasyCurve™ mop is an ideal choice for cleaning ceilings, walls and floors.

Features	Benefits
Hydrophilic polyurethane foam mop head	<ul style="list-style-type: none"> • Very absorbent and drip-free • Compatible with most common solutions, solvents and disinfectants
Laminated fabric face	<ul style="list-style-type: none"> • Clean and durable, for most surfaces
Ribbed microfiber face	<ul style="list-style-type: none"> • Better contaminant removal
Disposable flat-form mop head	<ul style="list-style-type: none"> • Economical and easy to install without tools • Will not pinch or tear gloves
Lightweight mop head and handle	<ul style="list-style-type: none"> • Weighs much less than roller-type and self-wringing mops • Significantly reduces operator fatigue
All components are autoclavable	<ul style="list-style-type: none"> • Ideal for a wide range of cleaning applications in the most critical sterile environments

Continued on reverse side



2743S

Features	Benefits
Curved, pivoting mop frame	<ul style="list-style-type: none"> • Provides effortless contact with surfaces • Ensures complete and uniform contact with all surfaces • Ensures proper “lift and pull” cleaning technique
Separate wringer in bucket, not on mop handle	<ul style="list-style-type: none"> • Removes weight from mop handle • Greatly reduces operator fatigue
Sieve wringer	<ul style="list-style-type: none"> • Easy and simple to use, no moving parts • Effectively wrings mop and removes contaminants from the mop head • Compatible with solutions, solvents, and disinfectants



Part No.	Description	Packaging
EasyCurve™ Mop and Hardware		
EZMP0300	EasyCurve disposable mop head, ribbed microfiber, 14 x 4.25 x 0.75" (36 x 11 x 2cm)	4/bag; 12 bags/case
EZMP0350	EasyCurve disposable mop head, ribbed microfiber, 14 x 4.25 x 0.75" (36 x 11 x 2cm), individually packaged, sterile	1/bag; 48 bags/case
2724	EasyCurve pivoting mop frame, stainless steel	1 each/case
2725E	QuickConnect interchangeable handle, 60" (150cm), stainless steel, lightweight, 1 piece, autoclavable	1 each/case
2643	QuickConnect interchangeable handle, 50" (130cm) extends to 92" (230cm), anodized aluminum, autoclavable	1 each/case
2645	QuickConnect interchangeable handle, 50" (130cm) extends to 92" (230cm), stainless steel, autoclavable	1 each/case
Bucket System		
2743S	Twin divided bucket w/stainless steel wringer; 5.3 gallon (20L) total volume, includes one cover, 4 non-autoclavable casters, polypropylene (seafoam)	1 each/case

Contec Healthcare offers many compatible mopping systems. Visit our website for details.