

Contec® Critical Site® Sterile Wipes

Sterile presaturated wipes designed for convenience in the pharmacy

Contec® Critical Site® Sterile Wipes are 100% polyester 4" x 4" (10.2 x 10.2 cm) low-lint wipes that contain 70% USP grade isopropyl alcohol and 30% USP grade purified water. The low-lint, nonwoven wipes are suitable for ISO Class 5 environments, compliant with Section 8.3 of USP Chapter <797>, and are exceptionally clean making them ideal for the most critical applications. Manufactured in Contec's ISO registered facility, the wipes are Validated Sterile per the Association for the Advancement of Medical Instrumentation (AAMI) 11137 Guidelines to a 10⁻⁶ Sterility Assurance Level (SAL). Sterility documentation is available upon request.

Presaturated Contec Critical Site Sterile Wipes create a repeatable, consistent technique that allows pharmacists and technicians to better control contamination during swabbing of critical entry points. These one-step wipes can be used in place of sterile individually wrapped alcohol prep wipes that require tearing of foil/paper pouches in the Primary Engineering Control (PEC) hoods. The low profile peel and reseal cleanroom-compatible packaging provides both quick access and easy closure, and requires minimal space in the PEC. All components of Contec Critical Site Wipes are made in the USA.



Contec Critical Site Sterile Wipes' 4" x 4" size make them ideal for use on smaller surfaces where precise application and sterility requirements are mandatory.

Features	Benefits
Presaturated Wipes	<ul style="list-style-type: none"> Contains consistent amount of 70% USP grade isopropyl alcohol (IPA) Allows for process repeatability Easy to grab wipe edges are ideal for dispensing
Validated Sterile	<ul style="list-style-type: none"> Features a Sterility Assurance Level of 10⁻⁶ per AAMI Guidelines Compliant with USP <797>
Synthetic, Low Lint Construction	<ul style="list-style-type: none"> Nonwoven, 100% low-lint polyester fabric is cleanroom processed Cleaner than sterile gauze and other similar products
Optimal Packaging	<ul style="list-style-type: none"> Peel and reseal packaging eliminates evaporation concern during storage when properly sealed Reduces potential risk of contamination from torn foil/paper 60 wipes per easily recognized pouch is ideal for pharmacy applications Minimal space requirements in the PEC. Little to no first air disruption due to the low profile packaging

Part no.	Description	Size	Packaging
HCPS6044	Contec Critical Site Sterile Wipes Polyester nonwoven wipes, 70% IPA	4" x 4" (10.2 x 10.2 cm)	60/pouch, 4 pouches/inner bag, 2 inner bags/outer bag, 3 outer bags/case; 1440 wipes/case