

## Wipe Guide

When it comes to choosing an effective wipe, Contec® Professional Lab Animal Science wants you to know you have options!

- Laundered rags often contain contamination that can interfere with your choice of chemistry
- Presaturated options can be expensive and tear easily, creating more contamination headaches.
- Paper and pulp-based products contain cellulose which can also interfere with chemical actives and reduce efficacy.

Premira® Microfiber Wipes, ContecClean™ Cloths, and Unifore® Surface Wipes are all made with 100% synthetic polymers, giving the broadest range of chemical compatibility.

Find the right wipe for the job and save money with increased efficiency, lowered cross contamination potential and reduced costs!



	Coverage	Cost	Cleaning Power
Premira® Microfiber Wipes	Best	Good	Best
ContecClean™ Cloths	Better	Better	Better
Unifore® Surface Wipes	Good	Best	Good

Durable all-around wipes with great coverage. Effective and efficient:

Part No.	Description	Size	Packaging	Saturation Unit	Recommended Sat Level (oz.)
PRMW0069	ContecClean Cloths, Perfed roll with 1 PRCN0001 canister	6" x 9" (15 x 23 cm) + 76 oz. (2.3 L) canister	60/roll; 22 rolls/case	roll	22
PRMW9069	ContecClean Cloths Perfed roll for use with PRCN0001 canister	6" x 9" (15 x 23 cm)	90/roll; 20 rolls/case	roll	32
PRMW0075	ContecClean Cloths Half fold, perfed roll for use with PRCN0007 bucket	10.5" x 12" (26 x 30.5 cm)	75/roll; 12 rolls/case	roll	60
PRMW1213	ContecClean Cloths Quarter folded for use with PRCN0005 bucket	12" x 13" (30.5 x 33 cm)	500/case	stack of 20	24
PRMW1217	ContecClean Cloths Quarter folded for use with PRCN0005 bucket	12" x 17" (30.5 x 43 cm)	150/box	stack of 20	32

Excellent mechanical cleaning power with patented microfiber. Picks up bedding, fur and debris with ease:

Part No.	Description	Size	Packaging	Saturation Unit	Recommended Sat Level (oz.)
PRMW3001	Premira Microfiber Wipes Quarter folded for use with PRCN0005 bucket	12" x 13" (30.5 x 33 cm)	600/case	stack of 25	32

Cost effective wipe ideal for applying your choice of chemistry:

Part No.	Description	Size	Packaging	Saturation Unit	Recommended Sat Level (oz.)
PRMW4008	Unifore Surface Wipes Perfed roll for use with PRCN0011 canister	10" x 12" (25.4 x 30.5 cm)	40/roll; 8 rolls/case	roll	18
PRMW4048	Unifore Surface Wipes Perfed roll for use with PRCN0011 canister	10" x 12" (25.4 x 30.5 cm)	40/roll; 48 rolls/case	roll	18
PRMW6030	Unifore Surface Wipes, perfed roll for use with PRCN0011 canister	10" x 12" (25.4 x 30.5 cm)	60/roll; 30 rolls/case	roll	22
PRMW1480	Unifore Surface Wipes, Quarter folded with PRCN0005 bucket and lid	12" x 13" (30.5 x 33 cm)	1,120/case	stack of 40	22
PRMW1680	Unifore Surface Wipes Quarter folded for use with PRCN0005 bucket	12" x 13" (30.5 x 33 cm)	1,280/case	stack of 40	22

For best results, it is recommended to replace wipes and chemistry after 7 days.