

## **EasyReach<sup>™</sup> Cleaning Tool**

One piece Stainless Steel EasyReach™ Cleaning Tool

## **Description**

The EasyReach Cleaning Tool is a one piece handle/ headframe used to gain access to hard to reach areas or confined spaces. The stainless steel construction is lightweight yet durable. Ideal to use inside isolators, fume hoods, glove boxes, RABS, etc. This unit is 304 stainless steel and can be autoclaved or cold sterilized. Autoclavable at 121°C for 30 minutes. Use with

Autoclavable at 121°C for 30 minutes. Use with all EasyReach cleaning pads: dry, presaturated, sterile, and non-sterile.

## **Technical Information**

Part No.	Description	Material
2665SF	Handle and frame, EasyReach, 17" long	Type 304 stainless steel



2665SF

## Notes:

- 1) Data shown are typical values and should not be used as product specifications.
- 2) The information presented here is applicable to the part numbers shown above, as well as to any product containing the same materials and produced under the same conditions, regardless of product size or packaging configuration. Please contact a Contec sales representative for more details.

www.contecinc.com

Revision date: 01/23/13

