



PROSAT[®] Wipes, Cellulose/polyester

Cellulose/polyester nonwoven wipe with 70% USP Grade IPA 30% DI Water

Description

These PROSAT wipes are made of strong cellulose/polyester nonwoven material that is very low in extractable residues. Saturated with a blend of 70% isopropyl alcohol and 30% DI water, these wipes are ideal for general purpose wiping and cleaning of equipment and work surfaces. The pouch opens easily to dispense wipes and seals after use to keep wipes saturated. Also available validated sterile for use in aseptic environments.

See part numbers on page 2

Technical Data

Wipe material	54% Cellulose / 46% Polyester	
Wipe construction	Hydroentangled staple fibers	
Solution of saturation; (blend ratio and constituents)	70% USP Grade Isopropanol, 30% Deionized Water	
Solvent Grade	USP	
Attribute; (units)	Value **	Test Method
Basis weight; (g/m ²)	53	Contec Method
Non-volatile residue, NVR		IEST-RP-CC004.3, Sec. 7.1.2
In deionized water; (g/m ²)	0.037	
In isopropanol; (g/m ²)	ND	
Specific ions		IEST-RP-CC004.3, Sec. 7.2.2
Sodium; (ppm)	11.9	
Chloride; (ppm)	18.2	
Particles, readily releasable		IEST-RP-CC004.2, Sec. 5.1
P ≥ 0.5µm; (x10 ⁶ /m ²)	18.6	
Fibers > 100µm; (x10 ³ /m ²)	21.7	IEST-RP-CC004.2, Sec. 5.2

** ND = None detected; levels are below detection limit of test equipment

Notes:

- 1) The information presented here is applicable to the part numbers shown above, as well as to any product containing the same materials and produced under the same conditions, regardless of product size or packaging configuration. Please contact a Contec sales representative for more details.
- 2) Data shown are typical values and should not be used as product specifications.
- 3) Valid product comparisons may only be obtained through side-by-side testing in the same test facility, under similar conditions.
- 4) Current and/or comparison data may be available. Please contact a Contec sales representative for details.

Revision date: 04/20/17
Page 1 of 2



Copyright © 2017 Contec, Inc. All rights reserved.

North America
P.O. Box 530
Spartanburg, SC 29304 USA
Phone: +1-864-503-8333
www.contecinc.com

Europe
Zl du Prat – RP 3707
56037 VANNES Cedex
FRANCE
+33 297 437 690

China
No.17 Longyun Road
Suzhou 215024
CHINA
+ 86-512-6274 4050

Part Number	Description	EA/PCH	PCH/OB1	OB1/CS	EA/CS
PSC20006	PROSAT wipes with 70% USP Grade IPA/30% DI Water, 9" x 11" (23 x 28cm)	50	4	7	1400

Part Number	Description	EA/PCH	PCH/OB1	OB1/OB2	OB2/CS	EA/CS
PSC20010	PROSAT wipes with 70% USP Grade IPA/30% DI Water, 9" x 9" (23 x 23cm)	50	7	1	3	1050

Part Number	Description	EA/PCH	PCH/OB1	OB1/OB2	OB2/CS	EA/CS
PSC20005	PROSAT wipes with 70% USP Grade IPA/30% DI Water, 9" x 11" (23 x 28cm), Sterile	50	7	1	4	1400
PSC20009	PROSAT wipes with 70% USP Grade IPA/30% DI Water, 9" x 9" (23 x 23cm), Sterile	30	9	1	3	810
PSCP0001	PROSAT wipes with 70% USP Grade IPA/30% DI Water, 6" x 9" (15 x 23cm), Sterile	50	4	2	5	2000

Key: EA = Each, PCH = Pouch, OB1 = Outer Bag 1, OB2 = Outer Bag 2, CS = Case

Notes:

- 1) The information presented here is applicable to the part numbers shown above, as well as to any product containing the same materials and produced under the same conditions, regardless of product size or packaging configuration. Please contact a Contec sales representative for more details.
- 2) Data shown are typical values and should not be used as product specifications.
- 3) Valid product comparisons may only be obtained through side-by-side testing in the same test facility, under similar conditions.
- 4) Current and/or comparison data may be available. Please contact a Contec sales representative for details.

Revision date: 04/20/17
Page 1 of 2



Copyright © 2017 Contec, Inc. All rights reserved.

North America
P.O. Box 530
Spartanburg, SC 29304 USA
Phone: +1-864-503-8333
www.contecinc.com

Europe
Zl du Prat – RP 3707
56037 VANNES Cedex
FRANCE
+33 297 437 690

China
No.17 Longyun Road
Suzhou 215024
CHINA
+ 86-512-6274 4050